

	No	Type		No	Type		No	Type
<b>AB Dental</b>			<b>Bionnovation</b>			<b>DEVA</b>		
Internal Hex 3,0	55021	11	Biomorse	55023	14	Bone Level 3.5	55064	56
Internal Hex 3.5 - 6.0	55049	41				Bone Level 4.5	55051	43
			<b>Biotech Dental</b>			<b>DIO</b>		
<b>ADIN</b>			Kontakt, Kontakt-S 3.6 - 4.8	55045	37	SM Submerged 3.8	55049	41
Touareg-S	55034	26				SM Submerged 5.0	55020	10
Touareg Closefit UNP 2.75	55059	51	<b>BlueSkyBio</b>			UF II Narrow	55064	56
			BIO Max NP Platform, All	55046	38	UF II Regular	55020	10
<b>AIRES</b>			BIO Max DP Platform 3.0-3.5	55046	38			
3.3	55049	41	BIO Max DP Platform 4.3-5.0	55037	29	<b>Easy Implant</b>		
			BIO Int Hex 3.5 Platform 3.7-4.1	55049	41	Ocean	55046	38
<b>Alpha-Bio</b>			BIO Int Hex 3.5 Platform 4.7	55051	43	Lilas	55037	29
Internal Hex	55034	26	BIO Int Hex 4.5 Platform 4.7-7.0	55051	43			
Conical Hex 3.2	55065	57	BIO One-stage Narrow, Regular, Wide	55014	4	<b>Edierre</b>		
Conical Hex 3.75-5.0	55051	43	BIO One-stage Short Regular & Wide	55014	4	3-lobe	55021	11
			BIO Quattro Narrow 3.3	55061	53	4-lobe	55025	16
<b>Alfa Gate</b>			BIO Quattro Regular 4.1 - 5.6	55062	54			
M+ NP 3.5	55044	36	BIO Conus 12, 3.0 Platform 3.0	55064	56	<b>Euroteknika</b>		
M+ RP 4.3 - 5.0	55037	29	BIO Conus 12, 3.5/4.0 Platform 3.5-4.0	55046	38	Internal Platform 3.0	55064	56
Porous, All diam	55034	26	BIO Conus 12, 4.5 x 6mm	55046	38	Internal Platform 3.5	55081	73
			BIO Conus 12, 4.5/5.0 Platform 4.5 - 5.0	55017	7			
<b>Alpha Dent</b>			BIO Conus 12, 5.0 x 6mm	55046	38	<b>Global D</b>		
IC, ICK	55067	59	BIO Trilobe 3.5	55016	6	In-Kone 4.0-5.0	55024	15
			BIO Trilobe 4.3-6.0	55017	7			
<b>Anker</b>			<b>Bredent</b>			<b>gmi</b>		
SB III 3.5	55023	14	SKY	55025	16	avantgard RP 3.75-4.25	55066	58
SB III 4.0-5.0	55020	10	copaSKY	55027	19	avantgard WP 4.75-5.75	55075	67
			mini2SKY	55033	25	frontier RP,WP	55051	43
<b>Ankylos</b>			fast & fixed Abutment	55083	75	phoenix NP 3.3	55038	30
C/X	55025	16				phoenix Std 3.3 - Wide 5.0	55011	1
			<b>BTI</b>			<b>Hahn</b>		
<b>Anthogyr</b>			Interna Universal	55051	43	Tapered 3.0	55058	50
Axiom Bone Level 3.4 - 5.2	55065	57	Interna Wide	55051	43	Tapered 3.5	55046	38
Axiom Tri-lobe	55065	57				Tapered 4.3	55046	38
Axiom Tissue Level Narrow, Regular	55050	42	<b>BTK</b>			Tapered 5.0	55037	29
			External Hex EN	55038	30	Tapered 7.0	55048	40
<b>Argon</b>			External Hex ER	55011	1	<b>Hiossen</b>		
K3pro Rapid 3.0 - 3.5	55026	17	External Hex EW	55011	1	HS System	55028	20
K3pro Rapid Platf 3mm	55019	9	Internal Hex IR	55049	41	ET System	55037	29
K3pro Rapid Platf 3mm Short	55044	36	Internal Hex IM	55051	43	ET 3.2	55047	39
			Internal Hex IW	55051	43	ET mini	55015	5
			Hexagonal Morse Taper KR	55081	73	<b>IDI Evolution</b>		
<b>Astra Tech</b>			Hexagonal Morse Taper KW	55056	48	Standard implant 3.75-4.0	55004	18
Osseospeed EV 3.0	55058	50	Octagonal Morse Taper SR	55014	4			
Osseospeed EV 3.6	55046	38	<b>Camlog</b>			<b>Implance</b>		
Osseospeed EV 4.2	55057	49	Screw-Line, Root-Line 3.3	55036	28	Tissue level	55014	4
Osseospeed EV 4.8	55055	47	Screw-Line, Root-Line 3.8, 4.3	55031	23	Bone level Regular	55037	29
Osseospeed EV 5.4	55048	40	Screw-Line, Root-Line 5.0, 6.0	55032	24	Bone level 3.3	55065	57
Osseospeed TX 3.0 S	55064	56	iSy	55065	57	<b>Implant Direct</b>		
Osseospeed TX 3.5-4.0 S	55046	38				Legacy 1/2/3/4 3.0/3.2	55012	2
Osseospeed TX 4.5-5.0 S	55017	7	<b>Ceralog</b>			Legacy 1/2/3/4 3.5/3.7 - 5.7/7.0	55051	43
			Hexalobe ceramic	55076	68	RePlant 3.5	55022	12
<b>Avinent</b>						RePlant 4.3 - 5.0	55010	13
External, Coral & Ocean 3.5	55038	30	<b>Conelog</b>			RePlant 6.0	55032	24
External, Coral & Ocean 4.1	55011	1	Screw-Line 3.3	55066	58	ReActive 3.5 Platform	55022	12
External, Coral 5.1	55013	3	Screw-Line 3.8, 4.3	55066	58	ReActive 4.3 - 5.0 Platform	55010	13
Internal, Coral & Ocean 3.5, 4.1	55013	43	Screw-Line 5.0	55075	67	Interactive 3.0 Platform	55046	38
Conical, Ocean 3.5	55066	58	<b>Dentack</b>			Interactive 3.4 Platform	55037	29
Conical, Ocean 4.1	55004	18	Internal Hex	55056	48	Scwindirect 5.0 PF 3.7 - 5.7	55070	62
						<b>Implant One</b>		
<b>BEGO</b>			3.3	55034	26			
Semados Mini	55012	2	3.75	55057	49	300 Series Std thread (blue)	55029	21
Semados S / RI 3.25-5.5	55034	26	4.5-5.5	55049	41	400 Series Std thread (purple)	55035	27
Semados RS / RSX 3.0	55082	74	3,8/4,5	55049	41	500 series Std thread (green)	55035	27
Semados RS / RSX 3.75-5.5	55034	26				300 Series Wide thread (blue)	55029	21
Semados SC / SCX 3.25-5.5	55034	26	<b>Dentegris</b>			400 Series Wide thread (purple)	55035	27
			3.3	55034	26	500 Series Wide thread (green)	55035	27
			3.75	55057	49			
			4.5-5.5	55049	41	<b>Implant Swiss</b>		
			3,8/4,5	55049	41	Tissue level	55014	4
<b>Bio3</b>						Bone level Regular	55037	29
Progressive 3.3, 3.8, 4.2, 5.0	55067	59	<b>C-Tech</b>			Bone level 3.3	55065	57
Advanced 3.3, 3.8, 4.2, 5.0	55067	59	EL	55067	59	<b>Implantvision</b>		
						Internal Hex 3.7-5.7	55034	26
<b>BioHorizons</b>			<b>Cowellmedi</b>					
Laser-Lok 3.0 Grey	55034	26	Inno Int Octa	55014	4			
Internal Laser-Lok 3.5 Yellow	55049	41	Inno Sub Hex	55037	29			
Internal Laser-Lok 4.5 Green	55051	43						
Internal Laser-Lok 5.7 Blue	55051	43	<b>C-Tech</b>					
External 4.0 Platform	55011	1	EL	55067	59			
<b>Biomet/3i</b>			<b>Dentium</b>					
External Hex	55011	1	Superline	55037	29			
Certain internal Purple	55036	28	NR Line	55066	58			
Certain internal Purple Platform 3.4	55066	58	Simple Line II GH 1,6	55014	4			
Certain internal Blue, Yellow, Green	55066	58						
T3 Purple	55036	28						
T3 Blue, Yellow, Green	55066	58						



---

	No	Type	No	Type
Zuga				
All	55027	19		